



ELK Biotechnology


LDHD Rabbit pAb

Catalog NO.: EA126

For research use only.

Overview

Product name	LDHD Rabbit polyclonal antibody
Source	Rabbit
Applications	WB
Species reactivity	Human, Rat, Mouse
Recommended dilutions	WesternBlot:1/1000-5000 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Recombinant Protein
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	55kDa
GeneID (Human)	197257
Human Swiss-Prot No.	Q85WU2
Cellular localization	Mitochondrial
Alternative Names	D Lactate Dehydrogenase antidoby, DLD antibody, Lactate dehydrogenase D antibody, MGC57726 antibody
Background	LDHD belongs to the D-isomer specific 2-hydroxyacid dehydrogenase family. D-lactate and D-glycerate dehydrogenase activities have been found in a similar protein in yeast.



Western blot analysis of 1) HeLa, 2) Mouse Skeletal Muscle , 3) Mouse Liver, 4) Rat Kidney with LDHD Rabbit pAb EA126 diluted at 1:5,000.