



ELK Biotechnology


SLUG Rabbit pAb

Catalog NO.: EA168


For research use only.

Overview

Product name	SLUG Rabbit polyclonal antibody
Source	Rabbit
Applications	WB, IHC
Species reactivity	Human, Rat, Mouse
Recommended dilutions	WesternBlot:1/500-1000 Immunohistochemistry:1/200-500 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Recombinant Protein
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	30kDa
GeneID (Human)	6591
Human Swiss-Prot No.	O43623
Cellular localization	Nucleus
Alternative Names	MGC10182 antibody, SLUGH1 antibody, SNAI2 antibody, Zinc finger protein SLUG antibody, WS2D antibody
Background	Transcriptional repressor. Involved in the generation and migration of neural crest cells.



Western blot analysis of 1) MCF7, 2) Mouse Heart Tissue, 3) Rat Heart Tissue, 4) Rat Brain Tissue with SLUG Rabbit pAb EA168 diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded human Breast carcinoma using SLUG (EA168) Rabbit pAb diluted at 1:500.