

EPG5 Rabbit pAb Catalog NO.: EA172 For research use only.

Overview

Product name EPG5 Rabbit polyclonal antibody

Source Rabbit

Applications WB, IHC

Species reactivity Human, Rat, Mouse

Recommended dilutions WesternBlot:1/1000-2000

Immunohistochemistry:1/200-500

NOTE: Optimal dilutions should be determined by the end user.

Immunogen Recombinant Protein

Species Human

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20° C. Avoid repeated freeze-thaw cycles.

Isotype IgG

Clonality Polyclonal

Concentration 1 mg/ml

Observed band 170-292kDa

GenelD (Human) 57724


Human Swiss-Prot No. Q9HCE0

Cellular localization cytoplasm


Alternative Names Ectopic P granules autophage protein 5 homolog antibody, HEEW1

antibody, hEPG5 antibody

Background N/A



Western blot analysis of Rat Liver Tissue with EPG5 Rabbit pAb EA172 diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Mouse Liver Tissue using EPG5(EA172) Rabbit pAb diluted at 1:500.