



ELK Biotechnology


TLR2 Rabbit pAb

Catalog NO.: EA184

For research use only.

Overview

Product name	TLR2 Rabbit polyclonal antibody
Source	Rabbit
Applications	IHC
Species reactivity	Human, Rat, Mouse
Recommended dilutions	Immunohistochemistry:1/200-500 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Recombinant Protein
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	N/A
GeneID (Human)	7097
Human Swiss-Prot No.	O60603
Cellular localization	Membrane
Alternative Names	CD282 antibody, TIL4 antibody, Toll like receptor antibody
Background	Toll-like receptor 2 also known as TLR2 is a protein that in humans is encoded by the TLR2 gene. TLR2 has also been designated as CD282. TLR2 is one of the toll-like receptors and plays a role in the immune system.



Immunohistochemical analysis of paraffin-embedded Mouse Spleen Tissue using TLR2 (EA184) Rabbit pAb diluted at 1:500.