



# ELK Biotechnology

TLR6 Rabbit pAb


Catalog NO.: EA185

For research use only.

## Overview

---

Product name	TLR6 Rabbit polyclonal antibody
Source	Rabbit
Applications	IHC
Species reactivity	Human, Rat, Mouse
Recommended dilutions	Immunohistochemistry:1/200-500 <b>NOTE: Optimal dilutions should be determined by the end user.</b>
Immunogen	Recombinant Protein
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	N/A
GeneID (Human)	10333
Human Swiss-Prot No.	Q9Y2C9
Cellular localization	Cell membrane, Cytoplasmic vesicle, Golgi apparatus
Alternative Names	CD286 antibody, Toll like receptor antibody
Background	Toll-like receptor 6 is a protein that in humans is encoded by the TLR6 gene. TLR6 has also been designated as CD286.



Immunohistochemical analysis of paraffin-embedded Mouse Spleen Tissue using TLR6 (EA185) Rabbit pAb diluted at 1:500.