

KCNN3(SK3) Rabbit pAb Catalog NO.: EA301 For research use only.

## Overview

Product name KCNN3(SK3) Rabbit polyclonal antibody

**Source** Rabbit

Applications IHC

Species reactivity Human, Rat

Recommended dilutions Immunohistochemistry:1/100-200

NOTE: Optimal dilutions should be determined by the end user.

Immunogen Synthetic Peptide

Species Human

**Storage** PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20° C. Avoid repeated freeze-thaw cycles.

**Isotype** IgG

**Clonality** Polyclonal

Concentration 1 mg/ml

Observed band 82kDa

GenelD (Human) 3782


Human Swiss-Prot No. Q9UGI6

Cellular localization Membrane

Alternative Names Kca2.3,hsk3,skca3

**Background** Forms a voltage-independent potassium channel activated by intracellular

calcium. Activation is followed by membrane hyperpolarization.



Immunohistochemical analysis of paraffin-embedded Rat BrainTissue using KCNN3(SK3) (EA301 ) Rabbit pAb diluted at 1:200.